

Direct Digital Pillow Speaker

- Modern slim case design
- ETL approved to UL1069
- Oxygen safe
- Easy to clean and disinfect
- Soft touch keypad
- Fluid resistant
- Audio jack included
- 7 vibrant colors now available

-Adding **COLOR** to the patient environment



NEW!

Quick Change Cable

1. Loosen screws



2. Remove cable




3. Install new cable



Pillow Speaker Ordering Guide

1. Determine Nurse Call System type i.e., Rauland, Hill-Rom, Jeron, etc.
2. Choose the appropriate connector type (listed below)



Dukane 18 Pin
ProCare



Dukane 18 Pin
Telligence




Executone 7 Pin
CareCom




Executone 16 Pin
Lifesaver



Fisher Berkeley 15 Pin
Lifesaver




Hill-Rom 20 Pin



Hill-Rom 37 Pin




Rauland 8 Pin
Responder 4




Rauland 9 Pin
Responder 5



Universal 1/4" Plug



West-Com® 10 Pin



Zettler® 15 Pin

3. Determine Make of Healthcare-Grade Television



LG

zenith



PHILIPS



SAMSUNG



TOSHIBA



4. Choose Pillow Speaker Case Color

